





Diagnostic reagent sample kits

As diagnostic tests move closer to the point-of-care, reagent storage and delivery methods become a challenge. LSO's East Coast location lab-on-chip reagent blisters allow reagents to be stored directly on the microfluidic cartridge, enabling the smaller footprint required by point-of-care diagnostic tests.



SAMPLE BLISTER FLUID VOLUME

DIMENSIONS	30µl	130µl	250µl	500µl	750µl	1,000µl
Height	0.108	0.175	0.190	0.26	0.245	0.268
Diameter	0.394	0.585	0.765	0.875	1.006	0.988
Flange	0.065	0.118	0.133	0.096	0.125	0.118
Pierce Diameter	0.264	0.375	0.500	0.625	0.750	0.752

Contact customer service at 1-603-692-9955 to order sample kits.

Credit Cards Accepted.

About LSO

J-Pac Medical is now part of Life Science Outsourcing. LSO is a leading FDA registered and ISO 13485 certified full-service medical device contract manufacturer. LSO offers clients services including cleanroom assembly, packaging, sterilization, medical device manufacturing, medical package testing, sterilization validation, fulfillment & distribution, cleaning & decontamination, and quality incubator services.



J-Pac Medical
Medical Device Outsourcing

Is now part of

Life Science Outsourcing

Over the past 25 years, Life Science Outsourcing has partnered with more than 1,400 medical device innovators to streamline production and navigate complex regulatory requirements.



Medical Device
Manufacturing



Medical Package
Testing



Sterilization &
Validation



Cleaning &
Decontamination



Fulfillment &
Distribution



Medical Device
Incubator Service